

# How to Download and Install Droplet Lab App for IOS?

App Installation · https://dropletlab.com/knowledge-doc/how-to-download-and-install-the-droplet-lab-mobile-application-for-ios/


1. Use the https://apps.apple.com/us/app/droplet-lab-sessile/id6445823509 and https://apps.apple.com/us/app/droplet-lab-pendant/id6445823559 URLs to land on the App Store screen with the application displayed.
2. To start download, tap the GET button. Doing so asks for the Face ID since it is the first time you are downloading the application.




3. Once verified, tap the OPEN button. This initiates the app installation.




4. Once installed, the app is displayed on the home screen. Tap the app icon to open the application.
5. As you tap the icon, the application asks for camera permissions. Tap the OK button.


"Sessile" Would Like to Access the Camera
Access to the camera is required to obtain image data for trials
Don't Allow
OK
Scan QR code to start using application